User Manual

Lumerus™ LED Exam Light

Strength in patient care.™

U.S. & Canada (888) 520-4998
International (801) 875-4998
www.mti.net

| | |
|--|----|
| APPLICABLE MODEL(S)..... | IV |
| REVISION HISTORY TABLE..... | IV |
| GENERAL INFORMATION..... | 1 |
| Purpose of this Manual..... | 1 |
| Intended Use..... | 1 |
| Symbols and Warnings..... | 2 |
| Symboles et Avertissements..... | 2 |
| Safety Instructions..... | 3 |
| Electromagnetic Interference..... | 3 |
| Electrical Requirements..... | 3 |
| Contraindications..... | 4 |
| Transportation and Storage Conditions..... | 4 |
| Disposal of Equipment..... | 4 |
| EQUIPMENT FEATURES..... | 4 |
| INSTALLATION..... | 6 |
| Required Tools..... | 6 |
| Unpacking..... | 6 |
| Installation..... | 7 |
| OPERATION..... | 8 |
| Operating Conditions..... | 8 |
| Controls..... | 9 |
| LIST OF ACCESSORIES..... | 10 |
| MAINTENANCE..... | 11 |
| Cleaning and Disinfection..... | 11 |
| Preventative Maintenance..... | 12 |
| Calls for Service..... | 12 |
| LIMITED WARRANTY..... | 12 |
| TROUBLESHOOTING GUIDE..... | 13 |
| SPECIFICATIONS LUMERUS™..... | 14 |

APPLICABLE MODEL(S)

This manual is for the following part numbers and serial numbers:

| PART NUMBER | DESCRIPTION | EFFECTIVE S/N |
|-------------|--|--------------------|
| 1000200 | Lumerus Exam Light w/Std 3-pronged Plug | 160310001- Current |
| 1000205 | Lumerus Exam Light w/IEC 320 Plug & Outlet Box | 160310001- Current |
| 1000207 | Lumerus Exam Light f/1" Posts | 160310001- Current |

REVISION HISTORY TABLE

| DATE | REVISION | CO#/ECO# | DETAILS |
|------------|----------|-----------------|--|
| 11/26/2014 | A | CO-0015/EC-0153 | Initial release |
| 4/7/2016 | B | CO-0115/EC-0359 | Updated instructions for improved neck design of lights. |
| 12/6/2019 | C | CO-0187 | Revised Proposition 65 Warning. Minor formatting fixes. Added sentence to Electromagnetic Interference section. Removed post mounted and caster base products. |

User Manual

MTI Lumerus™ LED Light

GENERAL INFORMATION

Purpose of this Manual

This manual covers the safety, installation, instructions for use, and maintenance of this equipment. It is intended that it be used by the service professional installing the equipment and by the medical professionals who have been trained to use it. No repair information is included with this manual because no operator repairs are intended.

Intended Use

This equipment is intended to be used as an examination light for general examinations and procedures performed by medical professionals.

PRODUCT INFORMATION

Date of Purchase: _____

Model Number: _____

Serial Number: _____

(see "EQUIPMENT FEATURES" on page 4 for location of serial number label)

Symbols and Warnings

This manual contains Notices, Cautions, and Warnings, which are meant to call attention to particular parts or sections of this manual. Since safety is our first concern, this manual should be read by all personnel who will be operating this medical equipment.

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WARNINGS/AVERTISSEMENTS

Are used to call attention to condition, practice, or procedure which could cause personal injury if not followed correctly. Please become familiar with all warnings and observe them at all times.

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Are used to call attention to a condition, practice, or procedure which could cause product damage if not used correctly.

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NOTICES/AVIS

Are used to call attention to helpful information regarding a condition, practice, or procedure.

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Safety Instructions

WARNING/AVERTISSEMENT

Read this manual before installation and use.

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WARNING/AVERTISSEMENT

No modification of this equipment is allowed.

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WARNING/AVERTISSEMENT

Do not operate equipment in the presence of flammable anesthetics.

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WARNING/AVERTISSEMENT

Assure controls, cords, and other parts of the equipment do not create a tripping or other safety hazard for operator or patient.

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WARNING/AVERTISSEMENT

This product can expose you to chemicals including Antimony Trioxide, which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov and see www.mti.net/faqs for details on California Proposition 65.

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Electromagnetic Interference

This equipment is designed and built to minimize electromagnetic interference with other devices. However, if interference is noticed between another device and this equipment, do one or more of the following:

1. Ensure other device is operating within its specifications.
2. Plug equipment into isolated circuit
3. Remove interfering device from room
4. Increase separation between equipment and interfering device
5. Contact MTI if electromagnetic interference persists

Electrical Requirements

WARNING/AVERTISSEMENT

To avoid risk of electrical shock, this equipment must be connected to a supply mains with protective earth.

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WARNING/AVERTISSEMENT

Grounding reliability can only be achieved when the equipment is connected to an equivalent receptacle marked "Hospital Only" or "Hospital Grade".

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WARNING/AVERTISSEMENT

Connecting electrical equipment to multiple socket-outlets effectively leads to creating a medical electronic system and can result in a reduced level of safety.

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WARNING/AVERTISSEMENT

Equipment must be positioned in normal use such that the power cords can be easily accessed to unplug.

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⚠ WARNING/AVERTISSEMENT

Equipment must be connected to an appropriate power source when loss of power source would result in an unacceptable risk.

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⚠ WARNING/AVERTISSEMENT

Do not use a two-three plug adapter to plug the equipment into a two conductor outlet, as grounding reliability is not maintained.

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Contraindications

None.

Transportation and Storage Conditions

| | | |
|---|----------------------|---|
|  | Ambient Temperature | +10°C to +40°C (+50°F to +104°F) |
|  | Relative Humidity | 30% to 75% (non-condensing) |
|  | Atmospheric Pressure | 700hPa to 1060hPa (0.69atm to 1.05atm) |


Disposal of Equipment

At the end of product life, the equipment, accessories, and other consumable goods may have become contaminated due to normal medical use. Consult local codes and ordinances for proper disposal of equipment, accessories and other consumable goods.


EQUIPMENT FEATURES**Description of Lumerus Exam Lights**

The Lumerus™ is an exam light that comes with a 3" light head. It can be mounted to a 1" post accessory or to the back of an MTI Exam Chair. The expected position of the operator is to the front of, to the side of, or behind the light. The expected position of the patient is below/underneath the light head.


LUMERUS™ LED LIGHT


1" POST ACCESSORY

LIGHT MOUNTING BRACKET MOUNTED TO LIGHT IN P/N 1000207

CHAIR MOUNTED LUMERUS™ LED LIGHT

INSTALLATION

CAUTION/PRÉCAUTION

To move equipment, lift only from the equipment's main structure.

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NOTICE/AVIS

If this product has been exposed to low temperatures, allow adequate time for it to warm to room temperature before putting it into service.

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NOTICE/AVIS

Inspect the shipping boxes for visible damage. Record on bill of lading, photograph, and report any damage to the carrier immediately. Request an inspection by the carrier if any concealed damage is noted.

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Required Tools

| | |
|---|--|
|  | Allen Wrenches, 3/32" (for 1" post mounted light) and 5/32" & 3/16" (for chair back mounted light) |
|  | Phillips Screw Driver (for chair back mounted light) |
|  | Utility Knife |


Unpacking

Open the boxes with the utility knife and remove all contents from the packaging.

Installation

The Lumerus™ can be purchased for installation on a 1" post accessory for your MTI Exam Chair or for installation on the back of an MTI Exam Chair.

1" Post Accessory Installation

1. Slide the light (with adapter) over the post and tighten the set screw enough to engage the groove, but not so tight that it prevents rotation of the light.


2. Plug the power cord into an appropriately marked receptacle.

! NOTICE/AVIS

Lumerus™ LED Exam lights can only mount to 1" posts that are designed to receive a light.

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To install the Lumerus™ with a standard 3-pronged plug to the back of an MTI Exam Chair

1. Remove the four bolts in the middle of the MTI Exam Chair back using a 5/32" Allen wrench.

2. Attach the Lumerus™ with the cable pointing away from the chair to the same four holes using the 4 screws provided with the Lumerus™ using a 3/16" Allen wrench.


3. Plug the power cord into an appropriately marked receptacle.

! WARNING/ AVERTISSEMENT

Disconnect the power to the MTI Exam Chair before installing the Lumerus™ LED Exam Light.

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To install the Lumerus™ with an IEC 320 plug to the back of an MTI Exam Chair


1. Remove the four bolts in the middle of the MTI Exam Chair back using a 5/32" Allen wrench.

2. Attach the Lumerus™ with the cable pointing away from the chair to the same four holes using the 4 bolts provided with the Lumerus™.

3. Remove the motor access panel under the seat of the MTI Exam Chair by removing two screws. Use a Phillips Screw Driver

4. Connect the 4-pin connectors of the wiring harness under the seat to the wiring harness of the electrical box assembly.

5. Install the electrical box assembly by attaching it with the same two screws previously removed.

6. Plug in the IEC plug of the Lumerus™ LED Light to the IEC receptacle on back of the electrical box. Push the plug in far enough to hear and feel a click to ensure full engagement of the plug

OPERATION

Operating Conditions

| | | |
|---|----------------------|--|
|  | Ambient Temperature | +10°C to +40°C (+50°F to +104°F) |
|  | Relative Humidity | 30% to 75% (non-condensing) |
|  | Atmospheric Pressure | 700hPa to 1060hPa (0.69atm to 1.05atm) |


Controls

⚠ CAUTION/PRÉCAUTION

Once the dimmer switch hits a stop in either the clockwise or counterclockwise direction, do not turn the dimmer switch further, otherwise the product may be damaged.

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On/Off and Brightness Control

1. To turn the light on, turn the ON/OFF switch clockwise. You will hear a click when the switch turns on.


2. To turn the light off, turn the ON/OFF Switch counter clockwise. You will hear a click when the switch turns off.

3. To brighten the light, turn the dimming switch clockwise to the desired brightness.

4. To dim the light, turn the dimming switch counterclockwise to the desired brightness.

Flexible Neck


The MTI Lumerus lights have a flexible neck to position the light head in the desired position.

1. To move the light:
A. Grasp the neck below the head of the light.
B. Bend the neck to the desired position.

Pivot Head

The MTI Lumerus lights have a pivot head, which has a 160° vertical and 180° horizontal rotation.

1. To pivot the light head:
A. Slightly loosen the pivot knob by turning it counterclockwise 1/8 of a turn.
B. Rotate the light head to the desired position.
C. Tighten the pivot knob by turning it clockwise.


⚠ CAUTION/PRÉCAUTION

Do not force the light head beyond the rotational stop or rotate the head more than 180° horizontally; otherwise, the light head will be damaged.

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Adjustable Bracket

The MTI Lumerus lights may have an adjustable bracket that allows a 90° pivot at the mounting location (1" post mounted and back mounted versions).

1. To adjust the bracket:
 - A. Loosen the knob on the bracket.
 - B. Move it to the desired position.
 - C. Tighten the knob.

LIST OF ACCESSORIES

| DESCRIPTION | USE | RESTRICTIONS |
|-------------|-----|--------------|
| None | | |

MAINTENANCE

WARNING/AVERTISSEMENT

There are no user serviceable parts inside. Refer all servicing to a qualified service technician.

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WARNING/AVERTISSEMENT

The Control of Hazardous Energy (lockout/tag out) §1910.147 procedure must be followed to disable equipment, preventing the release of hazardous energy while performing servicing and maintenance activities.

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CAUTION/PRÉCAUTION

This equipment may be damaged by solvents and dyes. Immediately remove any fluids spilled on the equipment.

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CAUTION/PRÉCAUTION

Any sprayed disinfectant or cleaning solution should be wiped up immediately after application and not let to stand. A small sample area should be tested before any large application of any cleaner/disinfectant.

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CAUTION/PRÉCAUTION

Avoid prolonged direct sunlight exposure to preclude color change of the plastic.

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NOTICE/AVIS

*MTI recommends cleaning the equipment frequently.
MTI recommande de nettoyer l'appareil fréquemment.*

Cleaning and Disinfection

MTI's standard policy on disinfecting MTI medical equipment is to use a 10% sodium hypochlorite (bleach) and 90% water solution. If stronger or other disinfecting and cleaning solutions are used, the surface is at risk to be damaged or possibly cause discoloring, cracking, or other types of damage to the surface. See current CDC Guideline for Disinfection & Sterilization in Healthcare Facilities.

MTI understands that many customers prefer to use a pre-mixed disinfecting solution in a spray or wipe. It is MTI's recommendation not to use sprays because over time the overspray will cause damage to the surfaces if not immediately wiped up. MTI recommends hand towels or wipes to avoid the possibility of overspray. Once the approved contact time has been obtained, remove and dry excess liquid remaining on the surface.

See "MTI Cleaning Solution Recommendation Letter" for more details on disinfecting solutions which have been tested on MTI's products. Available for download at <http://www.mti.net/downloads.aspx> in the Product Support Bulletins folder, found under the Resources tab on our website. Or, in the U.S. & Canada, call (888) 520-4998, or International, call (801) 875-4998 to request a copy.

Painted and Plastic Surfaces

All painted metal parts have been electrostatically painted and the paint is baked on for durability. These surfaces should be wiped down with a quality grade cleaner and a soft cloth or a damp cloth with mild soap once a week.

Preventative Maintenance

Periodically inspect the following areas:

1. All fasteners and pins should be in place and tightened securely.
2. Lubricate pins and moving joints as needed with lithium grease.
3. Power cord should be free of cuts or other visible damage.
4. All mechanical functions should operate properly.

Have an authorized service technician inspect your equipment once a year.

Calls for Service

Equipment's model and serial number information are required when calling for service.

Serial number is indicated by the  symbol located on the serial number tag.

If service is required, contact MTI

U.S. & Canada 888-520-4998

International 801-875-4998

Fax (801) 952-0548

customer.service@mti.net

Monday-Friday 8:00 am - 5:00 pm (MST)

www.mti.net

LIMITED WARRANTY

Medical Technology Industries, Inc. (hereinafter referred to as MTI), shall repair or replace products of its manufacture, which prove to be defective in material and/or workmanship for a period of twelve (12) months from the date of shipment to the customer (the "Warranty Period"). MTI will, at its option, provide parts to the customer or repair the defective part at MTI's factory or authorized repair facility. Except as it may otherwise specifically agree in writing, MTI shall not be liable for transportation or labor charges for repairs or adjustments or other work done on any of its products by MTI's authorized dealers or service organizations.

Customer shall notify MTI of any alleged claim or defect within fourteen (14) days from the date customer discovers, or upon reasonable inspection should have discovered, such alleged claim or defect (but in any event before the expiration of the Warranty Period). All repair parts are warranted for Ninety (90) calendar days from the date of shipment to the customer. Warranty is non-cumulative.

This warranty does not apply to any part or product which upon examination by MTI or its authorized agents indicates the product has been repaired or altered in any way by persons other than MTI or its authorized agents, has been misused, abused or used in a manner contrary to any instructions issued by MTI. Any parts or products claimed under this warranty must be properly packaged by customer and returned to MTI's factory with no liability to MTI for the parts, products or transportation charges thereon.

In no event, shall MTI have any liability for consequential, indirect, incidental, special, exemplary, punitive or enhanced damages, lost profits or diminution in value arising out of or relating to a product defect or the repair or replacement of any products manufactured or sold by it.

THE FOREGOING LIMITED WARRANTY IS MADE IN LIEU OF ALL WARRANTIES EXPRESSED OR IMPLIED, AND MTI MAKES NO OTHER WARRANTY WHATSOEVER WITH RESPECT TO THE PRODUCTS INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OR MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

TROUBLESHOOTING GUIDE

NOTICE/AVIS

For more information on this equipment, including circuit diagrams, parts and service manuals, and other descriptions, please contact MTI.

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| # | SYMPTOM | PROBLEM CAUSE | SOLUTION |
|---|-----------------------------|---|--|
| 1 | Light does not turn on. | Light is not plugged in. | Verify power cord is securely connected to facility power supply or to the IEC outlet on the back of the chair as applicable. |
| 2 | Light does not turn on. | Facility supply voltage | Check facility's power receptacle with a known working device. |
| 3 | Light does not turn on. | Facility supply voltage | Reset facility circuit breaker |
| 4 | Light does not turn on. | Electrical box fuses blown (Lumerus light with IEC plug only) | Check the fuses in the electrical box assembly. Remove by using a flat blade screw driver, press inward and turn counterclockwise. Only replace with fuses specified on electrical box assembly. |
| 5 | Lights do not dim/brighten. | Faulty dimmer switch | Contact Customer Service |

SPECIFICATIONS LUMERUS™

| DESCRIPTION | SPECIFICATION |
|---|---|
| Number of LEDs | 1 White LED |
| LED Life | 70,000 hours |
| Color Temperature (CCT) | 4,100° K , ±300° K |
| Working Distance | 16" (40.6 cm) |
| Light Intensity | 24,800 lux (2,300 fc) |
| Spot Size | 5" (12.7 cm) dia |
| Color Rendering Index (CRI) | 80 |
| Dimming Control | Infinitely adjustable 25-100% of full intensity |
| Electrical Input/Consumption | 100-240 VAC, 50-60Hz, 40W @ 3.3A |
| Head Rotation | 160° vertical and 180° horizontal |
| Mounting Options | Chair Back or 1" Post mounts |
| Power Cord 1" Post version Chair Back Mounted version | 3-wire grounded NEMA 5-15P Hospital Grade , 8' (2.44 m) long IEC 60320 C14 plug |
| Physical Dimensions: Head Diameter Length 32" Lumerus LED Length 42" Lumerus LED | 2.25" (5.7 cm) 40.6" (103 cm) Does not include mounting bracket. 52.8" (134.1 cm) Does not include mounting bracket. |
| Shipping Dimensions: Light - 32" Lumerus LED, IEC 320 Plug Light - 32" Lumerus LED, Std Plug Light - 42" Lumerus LED, Std Plug | 35" x 8" x 7" (88.9 cm x 20.3 cm x 17.8 cm) 5" x 24" x 25" (12.7 cm x 61 cm x 63.5 cm) 5" x 24" x 25" (12.7 cm x 61 cm x 63.5 cm) |
| Shipping Weight: Light (neck and head only) | 7.7 lbs (3.5 kg) |
| Warranty | 1-year limited |
| Listings | UL 153 UL 8750 CAN/CSA C22.2 No. 12 |

Strength in patient care.™

3655 W Ninigret Drive
Salt Lake City, Utah, 84104-6572
(800) 924-4655 • (801) 875-4999 • Fax (801) 952-0548
www.mti.net

© MTI, Inc. All rights reserved.
MTI reserves the right to make any
product changes without notice.

073092 Rev C