

WRAPPING

SUPPORT

# MTI Leg Wrapping Support\*

## A Breakthrough in Preventing Healthcare Professional Injury



- Allows caregiver to sit or stand
- Eliminates the difficulty of supporting patient's leg
- Maximizes patient comfort
- Allows one person to wrap most patient leg sizes
- Saves colleagues' time with healthcare provider cost savings
- Provides an ergonomic solution to potential caregiver injury

\*Patent Pending

# AVOID HEALTHCARE PROFESSIONAL INJURY

## MTI LEG WRAPPING SUPPORT\* CHAIR – BED – MOBILE

MTI's patent pending leg wrapping supports help reduce the amount of strain for both the patient and the health caregiver by supporting a patient's leg during the entire leg wrapping process.

This dual function support provides a calf support for wrapping the foot/heel and a heel support for wrapping the rest of the leg. The leg wrapping support is made from aluminum for a light-weight yet very strong support and is non-porous coated for ease in cleaning and disinfecting.



Foot/Heel wrapping position


Leg wrapping position

### MTI Chair-Mounted Leg Wrapping Support\*

The MTI Chair-Mounted Leg Support is attached to either side of the MTI 527, 527W, and 526 Series foot sections and can easily be moved to either side of the foot section to accommodate either the patient's right or left leg to easily wrap a leg or foot.

### Hospital Bed-Mounted Leg Wrapping Support\*

The Bed-Mounted Leg Wrapping Support is easily moved from bed to bed and is light weight making it portable throughout the entire facility.



Bed-mounted leg wrapping support, leg wrapping position shown

### Mobile Leg Wrapping Support\*

Light-weight yet strong, this portable support simplifies the wrapping of the leg or foot no matter where you use it. The mobile leg support folds-up into a compact unit for mobile applications. Take it in its carry bag with you anywhere: home, hospital, or on the road.



Foot/Heel wrapping position

Leg wrapping position

View Product  
Demonstration Video



Scan QR code with your smart phone to see a video demonstration of the leg wrapping support or go to:

[www.mti.net/MTI-Leg-Wrapping-Support.aspx](http://www.mti.net/MTI-Leg-Wrapping-Support.aspx)

For information  
or inquiries:

Call (800) 924-4655

Email [sales@mti.net](mailto:sales@mti.net)

Visit our web site at  
[www.mti.net](http://www.mti.net)

\*Patent Pending

© MTI, Inc. All rights reserved.  
MTI reserves the right to make any  
product changes without notice.

3655 W Ninigret Drive  
Salt Lake City, Utah 84104-6572  
800.924.4655 • 801.875.4999  
Fax 801.952.0548  
[www.mti.net](http://www.mti.net)

071051 Rev A (02/14)

Strength in patient care.™

