



G2

UPHOLSTERY COLOR GUIDE



Premium Upholstery

A beacon of strength, Premium Upholstery is both aesthetically pleasing and delightfully durable. Designed to complement any space and bring comfort to anyone in any MTI chair, stool, and side chair.


MTI's Premium Upholstery is FR free (fire-retardant free), compliant to California AB-2998 and has been formulated to be easily disinfected with sanitizing products containing sodium hypochlorite (common household bleach) diluted 1:10.




Taupe
215



Gunmetal
218



Mocha
237




Burgundy
238



Blue Ridge
221




Midnight Blue
239



Black
217



Slate
250



Dove
241



Premium
Touch Sample

Blue Ridge
221

Ultraleather[®] Upholstery

A premier in the world of upholstery fabrics and unmatched in comfort and high-wear durability. Thoughtfully developed with high-strength fibers to create a leading earth-friendly material free from PVC, heavy metals, stabilizers, VOCs, POPs, and Formaldehyde. The soft and supple hand feel is designed to delight while defining the ultimate in luxury in a human and eco-friendly way.

Ultraleather[™]
Touch Sample



Arctic
242



Papyrus
223



Pine Cone
243



Papyrus
223



Garnet
244



Diplomat Blue
236



Raven Wing
225



Schooner
226



Dove Grey
245



Silver Pearl
246

Specifications¹

	Premium Upholstery	Ultraleather [®] Upholstery
California Technical Bulletin 117-2013	•	•
PVC Free		•

Cleaning & Disinfection

MTI's standard policy on disinfecting medical equipment is to use a solution with a 10% sodium hypochlorite (bleach) and 90% water (1:10 ratio)². For premixed disinfectant options, see MTI Suggested Disinfectants Policy downloadable at www.mti.net/DisinfectantsPolicy

¹ It is the customer's responsibility to order products to meet local regulations

² See current CDC guidelines for Disinfection and Sterilization in Healthcare Facilities

Colors in this brochure are as closely represented as possible. However, due to the nature of monitors, printed material, and dye lots, the nature of the actual products may differ slightly from how they appear.

For more information contact MTI. Call us at (800) 924-4655 (US & Canada) or (801) 875-4999 (International) or visit us at www.mti.net.

MTI, Inc.
 3655 W. Ninigret Drive • Salt Lake City, Utah, 84104-6572
 Phone: 800-924-4655 • 801-875-4999 • Fax 801-952-0548
www.mti.net

Specialized Equipment
 for Specialty Healthcare
Medical Technology Industries

